

Product Name:	ABS-DDL2-18
	<p>           USB Digital Data Logger (DDL), complied to the current CDC vaccine storage monitor guidelines            2 years certification of calibration            “buffered” probe in the product simulated solution            °F/°C Switchable Readout            Min/Max historical temp point            Detachable Dual Probe            High and Low Temp Alarm            Visual Alarm            Audible Alarm            Operation Temp Range: -50.00 to 60.00°C (-58 to 140°F)            Accuracy of +/-0 .25° C            Readout resolution: 0.01            Low-Battery Indicator            Temp Data downloaded via USB drive (not included)            User programmable Logging Interval (recommend no less than every 30 mins)            Store up to 3 years of data (30-minute intervals)         </p>
Description:	