



Sliding and Fixed Baskets
Freestanding 10.5 and 5.2 ft³ Refrigerators

ABT-SLD-BSKT-FULL-5.2




Compatible Refrigerator Size	5.2 ft ³
Max Quantity	3
Basket type	Sliding
Width x Depth x Height	16 7/8" x 13 5/8" x 3"
Max Load	30 lbs.

ABT-BSKT-BASE-FIXED-5.2




Compatible Refrigerator Size	5.2 ft ³
Max Quantity	1
Basket type	Fixed
Width x Depth x Height	18 3/8" x 8 1/2" x 6 1/4"
Max Load	44 lbs.

ABT-SLD-BSKT-FULL-10.5



Compatible Refrigerator Size	10.5 ft ³
Max Quantity	6
Basket type	Sliding
Width x Depth x Height	16 7/8" x 15 3/4" x 3"
Max Load	35 lbs.

ABT-BSKT-BASE-FIXED-10.5



Compatible Refrigerator Size	10.5 ft ³
Max Quantity	1
Basket type	Fixed
Width x Depth x Height	18 3/8" x 10 5/8" x 6 1/4"
Max Load	55 lbs.

