



ABT-DWRKIT-PL23*
ABT-DWRKIT-PL49*

Aerial View



Front View



Compatible Models

- ABT-HC-PL-23
- ABT-HC-SPL-23
- ABT-HC-PL-49
- ABT-HC-SPL-49
- ABT-HC-PLF-23
- ABT-HC-PLF-49

Drawers Per Door	5
Drawer Type	Sliding
Width x Depth x Height	18 3/4" x 23 3/4" x 2"
Max Load	65 lbs.

*Factory Install only

