

Product Name:	ABS-WF-DDL-18
	<p>WiFi Digital Data Logger (DDL), complied to the current CDC vaccine storage monitor guidelines</p> <p>2 years certification of calibration</p> <p>"buffered" probe in the product simulated solution</p> <p>°F/°C Switchable Readout</p> <p>Min/Max historical temp point</p> <p>Detachable Probes (2 included. Can be used with only one probe)</p> <p>High and Low Temp Alarm</p> <p>Visual Alarm</p> <p>Audible Alarm</p> <p>Operation Temp Range: -50.00 to 60.00°C (-58 to 140°F)</p> <p>Accuracy of +/-0.25° C</p> <p>Readout resolution: 0.01</p> <p>Low-Battery Indicator</p> <p>Temp Data downloaded via WiFi network. Easy to set up.</p> <p>User programmable Logging Interval (recommend no less than every 30 mins)</p> <p>Store up to 3 years of data (30-minute intervals)</p>
Description:	