

Product Name:	ABS-WF-DDL-18
	<p>WiFi Digital Data Logger (DDL), complied to the current CDC vaccine storage monitor guidelines 2 years certification of calibration “buffered” probe in the product simulated solution °F/°C Switchable Readout Min/Max historical temp point Detachable Probes (2 included. Can be used with only one probe) High and Low Temp Alarm Visual Alarm Audible Alarm Operation Temp Range: -50.00 to 60.00°C (-58 to 140°F) Accuracy of +/-0 .25° C Readout resolution: 0.01 Low-Battery Indicator Temp Data downloaded via WiFi network. Easy to set up. User programmable Logging Interval (recommend no less than every 30 mins) Store up to 3 years of data (30-minute intervals) Connects to 2.4GHz Radio Frequency</p>
Description:	