

Product Name:	ABT-Wifi1-DL
Description:	<p>Single (or Dual) Probe WiFi Data Loggers are high speed wireless modules with PEAPv0 enterprise security, capable of collecting, storing and transmitting data wirelessly, connects to 2.4GHz frequency. The transmitter passes information to a standard access point which can be accessed by an Wi-Fi enabled network. Each transmitter monitors against preset conditions that are defined by the user and can provide audio and visual alerts. Additional alerts can be provided through a variety of methods such as SMS, text messages, voice, pager, cell phone, fax and email. Information recorded in °F or °C to the database is time-stamped and cannot be altered through the user interface.</p> <ul style="list-style-type: none"> • Alarms: Visual and audio alarm indicators – can be silenced locally • Housing: Easy zip tie mount hooks with screw holster wall mount for each mount and removal • Temperature Range: -4 °F to 194 °F (-20 °C to 90 °C) • Tolerance: Digital Temperature Probe ±0.5 °C • Ambient Operating Range: 32 °F to 167 °F (0 °C to 75 °C). Up to 95% RH • Power Supply: (4) 1.5V AA Lithium Batteries (included). Any AA 1.5V batteries can be used (lithium or alkaline) • Probe: Cable Length - 12 Ft.; Probe Length – 3.6" • Connects to 2.4GHz Radio Frequency <p>Product Dimensions: 4.25" H x 3.88"W x 1.38" (thickness) Shipping Dimensions: 8"x 8"x 8"; 2lbs.</p>