ABT-HC-CS-16
Standard Glass Door Chromatography Refrigerator 16 Cu. Ft. capacity One swing glass door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Three adjustable shelves (adjustable in ½" increments) Forced draft air circulation Heavy duty door handle and keyed door lock One GFCI protected duplex electrical outlet One 2" access port with spring loaded, self-closing gasketed metal cover Magnetic door gasket for positive seal LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (¾") Advanced hydrocarbon detector White powder-coated steel exterior and interior High density urethane foam cabinet insulation Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 25 3/4" W x 29 ¾" D x 79" H with casters; 75 ½" H without casters Cycle defrost 115V, 60 Hz, 3 Amps, 1/5 HP ENERGY STAR CERTIFIED ETL listed Shipping weight: 283 lbs. Note: This unit requires two (2) outlets