Product Name:	ABT-HC-CS-23
Description:	Standard Glass Door Chromatography Refrigerator 23 Cu. Ft. capacity One swing glass door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door lock One GFCI protected duplex electrical outlet One 2" access port with spring-loaded, self-closing gasketed metal cover Magnetic door gasket for positive seal Heavy duty door handle Four adjustable shelves (adjustable in ½" increments) Forced draft air circulation LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (¾") Advanced hydrocarbon detector White powder-coated steel exterior and interior High density urethane foam cabinet insulation Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 27 5/8" W x 34 ¾" D x 81 ¾" H with casters; 78 ¼" H without casters Cycle defrost 115V, 60 Hz, 3 Amps, 1/5 HP ENERGY STAR CERTIFIED ETL listed
	115V, 60 Hz, 3 Amps, 1/5 HP ENERGY STAR CERTIFIED