

Product Name:	ABT-HC-CS-49
Description:	<p>Standard Glass Door Chromatography Refrigerator 49 Cu. Ft. capacity Two swing glass doors, self-closing, right hinged Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door locks One GFCI protected duplex electrical outlet One 2" access port with spring-loaded, self-closing gasketed metal cover Heavy duty door handles Magnetic door gasket for positive seal Eight adjustable shelves (adjustable in ½" increments) Forced draft air circulation LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (3/8") Advanced hydrocarbon detector White powder-coated steel exterior and interior High density urethane foam cabinet insulation Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 54 3/4" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters Cycle defrost 115V, 60 Hz, 4.5 Amps, 1/4 HP ENERGY STAR CERTIFIED ETL listed Shipping weight: 538 lbs. Note: This unit requires two (2) outlets</p>