Product Name:	ABT-HC-CS-69
Description:	Standard Glass Door Chromatography Refrigerator 69 Cu. Ft. capacity Three sliding glass doors, self-closing (energy efficient, double pane vacuum insulated) Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door locks One GFCI protected duplex electrical outlet One 2" access port with spring-loaded, self-closing gasketed metal cover Integrated door handles and magnetic door gasket for positive seal Twelve adjustable shelves (adjustable in ½" increments) Forced draft air circulation LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (38") Advanced hydrocarbon detector White powder-coated steel exterior and interior High density urethane foam cabinet insulation Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 81 3/4" W x 33 1/4" D x 78 ¼" H without casters; 81 ¾" H with casters Cycle defrost 115V, 60 Hz, 7.5 Amps, 1/2 HP ETL listed Shipping weight: 738 lbs. Note: This unit requires two (2) outlets