Product Name:	ABT-HC-CP-16-TS
Description:	TempLog Premier Glass Door Chromatography Refrigerator 16 Cu. Ft. capacity
	One swing glass door, self-closing, right hinged, not reversible Microprocessor temperature controller
	Adjustable operating temperature range: 1°C - 10°C
	Touch Screen Digital Temperature Display and Alarm Module
	Password-protected control panel
	Sample simulated digital temperature display; .1° Resolution, adjustable
	temperature display in Celsius or Fahrenheit
	Active temperature scale bar display Resettable 1-point min/max history and digital 24-hour temperature chart
	Temperature data logging
	USB port for data transfer and software updates
	Real-time operational display clock
	4-20mA output
	Electronic troubleshooting guide
	Electronic QR Code to access Owner's Manual & technical support contact info
	Adjustable alarm volume
	Audible and visual high/low temperature alarms 2 probes; 1 in air and 1 in sample bottle
	Sensor error alarm, power failure alarm, door ajar alarm, low battery alarm
	Operation indicator
	Alarm validation (audible, visual, and remote alarm) and alarm mute with ring
	back Remote alarm contacts Three adjustable shelves (adjustable in 1/" ingrements)
	Three adjustable shelves (adjustable in ½" increments) Forced draft air circulation
	Keyed door lock
	One GFCI protected duplex electrical outlet
	Two 2" access ports with spring-loaded, self-closing gasketed metal covers
	Magnetic door gasket for positive seal
	Heavy duty door handle
	LED interior lights are safety shielded and switch controlled
	Swivel Casters (3 ½") Probe access port (3/4")
	Advanced hydrocarbon detector
	White powder-coated steel exterior and interior
	High density urethane foam cabinet insulation
	Hydrocarbon, natural refrigerant (R290)
	Three year parts and labor warranty, plus an additional two year compressor parts warranty
	Exterior dimensions: 26 1/2" W* x 29 3/4" D x 79" H with casters; 75 ½" H without casters. *Width w/access ports closed; 34 1/2" width w/access ports open
	Cycle defrost
	115V, 60 Hz, 3 Amps, 1/5 HP
	ENERGY STAR CERTIFIED
	ETL listed
	Shipping weight: 283 lbs. Note: This unit requires two (2) outlets