Product Name:	ABT-HC-CP-49-TS
Description:	Templog Premier Glass Door Chromatography Refrigerator 49 Cu. Ft. capacity
	Two swing glass doors, self-closing, right hinged, not reversible Microprocessor temperature controller
	Adjustable operating temperature range: 1°C - 10°C
	Touch Screen Digital Temperature Display and Alarm Module
	Password-protected control panel Sample simulated digital temperature display; .1° Resolution, adjustable temperature
	display in Celsius or Fahrenheit
	Active temperature scale bar display
	Resettable 1-point min/max history
	Digital 24-hour temperature chart
	Temperature data logging USB port for data transfer and software updates
	Real-time operational display clock
	4-20mA output
	Electronic troubleshooting guide
	Electronic QR Code to access Owner's Manual & technical support
	contact information
	Adjustable alarm volume
	Audible and visual high/low temperature alarms
	2 probes; 1 in air and 1 in sample bottle
	Sensor error alarm, power failure alarm, door ajar alarm, low battery alarm
	Operation indicator
	Alarm validation (audible, visual, and remote alarm) Alarm mute with ring back
	Remote alarm contacts
	Eight adjustable shelves (adjustable in $\frac{1}{2}$ " increments)
	Forced draft air circulation
	Keyed door locks
	One GFCI protected duplex electrical outlet
	Two 2" access ports with spring-loaded, self-closing gasketed metal covers Magnetic door gasket for positive seal and heavy duty door handles
	LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½")
	Probe access port (3/4")
	Advanced hydrocarbon detector
	White powder-coated steel exterior and interior
	High density urethane foam cabinet insulation
	Hydrocarbon, natural refrigerant (R290)
	Three year parts and labor warranty, plus an additional two year compressor parts warranty
	Exterior dimensions: 55 1/2" W* x 34 3/4" D x 81 ³ / ₄ " H with casters; 78 ¹ / ₄ " H without casters. *Width w/access ports closed; 63 1/2" width w/access ports open
	Cycle defrost
	115V, 60 Hz, 4.5 Amps, 1/4 HP
	ENERGY STAR CERTIFIED ETL listed
	Shipping weight: 538 lbs. Note: This unit requires two (2) outlets