

Product Name:	ABT-HC-UCBI-0420
Description:	<p>Premier Undercounter Freezer Built-In 4.2 Cu. Ft. Capacity One solid swing door, right hinged* Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in glycol bottle Adjustable operating temperature range: -15°C – -25°C Two fixed shelves Keyed door lock Recessed door handle Magnetic door gasket for positive seal Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 23 ¾"W x 24 ½"D x 33 3/8"H Manual defrost 115V, 60 Hz, 1.5 Amps 1/8 HP Energy Consumption 0.97 kWh/day UL/C-UL listed ENERGY STAR CERTIFIED Shipping weight: 132 lbs. *This product is available left hinged, product number ABT-HC-UCBI-0420-LH</p>