ABT-HC-UCBI-0404-ADA-LH
Premier Undercounter Refrigerator Built-In  1.6 Cu. Ft. Capacity  One solid swing door, left hinged  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  P. probes; 1 in air and 1 in glycol bottle  Adjustable operating temperature range: 1°C – 10°C  Three shelves (two adjustable/one fixed) with guard rail on back to  prevent product from pushing against back wall  Back wall evaporator cover  Forced draft circulation  Keyed door lock  Magnetic door gasket for positive seal  Heavy duty door handle  Leveling legs  LED interior light  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Two year parts and labor warranty, plus an additional three year  compressor parts only warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 23 ¾ "W x 24 ½" D x 31 15/16 "H  Dycle defrost  15V, 60 Hz, 1.74 Amps, 1/8 HP  ENERGY STAR CERTIFIED  ADA compliant  JL/C-UL listed  Shipping weight: 150 lbs.