

Product Name:	ABT-HC-UCFS-0104
Description:	<p>Premier Countertop Freestanding Refrigerator  1 Cu. Ft. Capacity  One solid swing door, right hinged*  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  2 probes; 1 in air and 1 in glycol bottle  Adjustable operating temperature range: 1°C – 10°C  2 adjustable shelves  Back wall evaporator cover  Forced draft circulation  Keyed door lock  Integrated door handle  Magnetic door gasket for positive seal  Leveling legs  LED interior light  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 17 ¼" W x 19 ¼" D x 21 ¼"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Energy Star Certified for Laboratory Refrigerators and Freezers  Cycle defrost  115V, 60 Hz, .90 Amps, 1/8  Energy Consumption 0.71 kWh/day  UL/C-UL listed  ENERGY STAR CERTIFIED  Shipping weight: 91 lbs.  *This product is available left hinged, product number ABT-HC-UCFS-0104-LH</p>