

Product Name:	ABT-HC-RFC-15A
Description:	<p>Standard Hydrocarbon Refrigerator & Freezer Combination ADA Compliant 14 cu. ft. capacity (10 ref./4 freezer) Two exterior solid swing doors, right hinged, reversible Mechanical thermostat Refrigerator temperature range: 0°C to 10°C Freezer chamber: -15°C to -10°C with two auto defrost cycles warming to 0°C per day Two adjustable glass shelves (refrigerator) and two-gallon door bins, plus one half gallon door bin 1 full width glass shelf (freezer) and 2 fixed door bins Magnetic door gasket for positive seal Integrated door handles Leveling legs Interior LED refrigerator light Two probe access ports (3/8"); one, refrigerator/one, freezer White textured cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 27 3/8"W x 29"D x 60 3/8"H Auto defrost 115V, 60 Hz, 1.1 Amps, 1/10 H.P. Shipping weight: 135 lbs. DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assist in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.</p>