

Product Name:	ABT-HC-FRFC10
Description:	<p>Premier Flammable Refrigerator & Freezer Combination 9 Cu. Ft. capacity (5 Ref./4 Freezer) Two exterior solid swing doors right hinged, not reversible Microprocessor temperature controller, refrigerator and freezer Digital temperature display, refrigerator and freezer Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C Audible and visual high and low temperature alarms, refrigerator and freezer 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer Remote alarm contacts, refrigerator and freezer 3 adjustable shelves refrigerator/2 fixed shelves (with 3 removable baskets) freezer Keyed door locks Magnetic door gasket for positive seal Recessed door handle, refrigerator and freezer Leveling legs Stacking kit Two probe access ports (3/8"); one, refrigerator/one, freezer White powder coated exterior HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty Exterior dimensions: 23 5/8"W x 24 1/4"D x 65 3/8"H Manual defrost freezer / cycle defrost refrigerator 115V, 60 Hz, 1.2 Amps Refrigerator and 1.4 Amps Freezer, 1/9 HP Refrigerator and 1/8 HP Freezer UL/C-UL listed Units meet NFPA guidelines 45 and 70 Shipping weight: 205 lbs. *Please note this unit requires 2 outlets.</p>