

| | |
|---------------|--|
| Product Name: | ABT-FRS-30 |
| Description: | <p>General Purpose Flammable Storage Refrigerator 30 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat Adjustable operating temperature range: 1°C – 10°C Five adjustable shelves Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle Leveling legs Probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging One year parts and labor warranty Exterior dimensions: 35 3/4"W x 35 1/2"D x 72 1/2"H Cycle defrost 115V, 60Hz, 6 Amps, 1/3 HP UL/C-UL listed Units meet NFPA guidelines 45 and 70 Shipping weight: 381 lbs.</p> |