Product Name:	ABT-EFS-30
Description:	General Purpose Hazardous Location (Explosion Proof) Freezer 30 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat sealed in UL listed explosion proof box Adjustable operating temperature range: -15°C – -25°C Four fixed shelves Keyed door lock Magnetic door gasket for positive seal Stainless steel door handle Leveling legs White textured cabinet and door High density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging Exterior dimensions: 35 ³/4"W x 35"D x 77 ¹/2"H Manual defrost Not equipped with line cord and plug; must be hardwired in conduit to 120 VAC power source 115V, 60Hz, 6 Amps, 1/3 HP UL/C-UL listed One year parts and labor warranty Meet NFPA and OSHA guidelines for hazardous locations, Class 1, Division II, Group C & D. These units are not intended for use in general clinical or research laboratories. Shipping weight: 390 lbs.