Co 30 Or no Me Ac Fr Th Ke	General Purpose Hazardous Location (Explosion Proof) Refrigerator/Freezer Combination 0 Cu. Ft. Capacity (25 Ref./5 Freezer) One solid swing door/one interior solid swing door-freezer chamber only, right hinged, ot reversible Mechanical thermostat sealed in UL listed explosion proof box Adjustable operating temperature range: Refrigerator: 1°C – 10°C Greezer: -15°C – -25°C Three adjustable refrigerator shelves Keyed door lock
Description: Description: Hi CF Ex Ma No to 11 UL Or Ma &	Magnetic door gasket for positive seal Stainless steel door handle eveling legs White textured cabinet and door ligh density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging Exterior dimensions: 35 <sup>3</sup> / <sub>4</sub> "W x 36 <sup>1</sup> / <sub>4</sub> "D x 77 <sup>1</sup> / <sub>2</sub> "H Manual defrost lot equipped with line cord and plug; must be hardwired in conduit to 120 VAC power source 15V, 60Hz, 6 Amps, 1/3 HP JL/C-UL listed One year parts and labor warranty Meet NFPA and OSHA guidelines for hazardous locations, Class 1, Division II, Group C a D.