Product Name:	ABT-HC-19-LH
Description:	
	19 Cu. Ft. capacity
	One swing glass door, self-closing (energy efficient, double pane
	vacuum insulated), LEFT hinged, not reversible
	Microprocessor Temperature Controller Adjustable operating temperature range: 1°-10°C
	Standard Digital Temperature Display and Alarm Module with
	battery back-up to support temperature monitoring and alarms
	Sample simulated digital temperature display; .1° Resolution
	Adjustable temperature display in Celsius or Fahrenheit
	Audible and visual high/low temperature alarms
	Sensor error alarm
	Power failure alarm
	Resettable 1-point Min/Max History
	Remote alarm contacts
	Alarm validation (audible, visual, and remote alarm)
	Alarm mute with ring back
	Four adjustable shelves (adjustable in ½" increments) Forced draft air circulation
	Keyed door lock
	Magnetic door gasket for positive seal
	Recessed door handle
	LED interior lights are safety shielded and switch controlled
	Swivel casters (adds 4" to height)
	Probe access port (3/8")
	White laminated vinyl exterior
	White aluminum liner and stainless steel floor interior
	HFC-free refrigerant (R290), foam insulation and packaging
	Two year parts and labor warranty, plus an additional three year
	compressor parts warranty
	Exterior dimensions: 27 1/8" W x 25 1/2" D x 82 5/8" H with casters; 78 5/8" H without casters
	Cycle defrost
	115V, 60 Hz, 5.4 Amps, 1/3 HP, BTU/HR 4932
	UL/C-UL listed
	Shipping weight: 295 lbs.