

Product Name:	ABT-HC-RFP-17
Description:	<p>Standard Natural Refrigerant Laboratory Refrigerator 17 Cu. Ft. capacity One solid swing door, right hinged, not reversible Microprocessor Temperature Controller Adjustable operating temperature range: 1°-10°C Standard Digital Temperature Display Sample simulated digital temperature display; .1° Resolution Audible and visual high/low temperature alarms Remote alarm contacts 5 adjustable shelves Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle Leveling legs One probe access port (3/4") White textured cabinet and door High density urethane foam cabinet and door insulation Natural, Hydrocarbon Refrigerant (R600a) One year parts and labor warranty, plus an additional four-year compressor parts warranty Exterior dimensions: 30 3/4" W x 28 5/8"D x 67 3/4"H Cycle defrost 115V, 60 Hz, 2.5 Amps, 1/4 HP ENERGY STAR CERTIFIED ETL listed Shipping weight: 240 lbs.</p>