Product Name:	ABT-HC-SLS-12
Product Name: Description:	ABT-HC-SLS-12 Standard Solid Door Laboratory Refrigerator 12 Cu. Ft. capacity One swing solid door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: $1^{\circ}C - 10^{\circ}C$ Keyed door lock Heavy duty door handle Magnetic door gasket for positive seal Three adjustable shelves (adjustable in $\frac{1}{2}$ " increments) Forced draft air circulation LED interior lights are safety shielded and switch controlled Swivel Casters ( $3\frac{1}{2}$ ") Probe access port ( $3/8$ ") White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 25" W x 29 $\frac{3}{4}$ " D x 65 $\frac{3}{4}$ " H with casters; 62 $\frac{1}{4}$ " H without casters Cycle defrost
	115V, 60 Hz, 3 Amps, 1/5 HP ETL listed Shipping weight: 227 lbs.