Product Name:	ABT-26S-TS
Description:	TempLog Premier Solid Door Laboratory Refrigerator
Description.	26 Cu. Ft. capacity
	One swing solid door, self-closing, right hinged, not reversible
	Microprocessor temperature controller
	Adjustable operating temperature range: 1°C - 10°C
	Touch Screen Digital Temperature Display and Alarm Module with battery back-up to
	support temperature monitoring and alarms
	Password-protected control panel
	Sample simulated digital temperature display; .1° Resolution, adjustable temperature
	display in Celsius or Fahrenheit
	Active temperature scale bar display
	Resettable 1-point min/max history Digital 24-hour temperature chart
	Temperature data logging
	USB port for data transfer and software updates
	Real-time operational display clock
	4-20mA output
	Electronic troubleshooting guide
	Electronic QR Code to access Owner's Manual & technical support
	contact information
	Adjustable alarm volume
	Audible and visual high/low temperature alarms
	2 probes; 1 in air and 1 in sample bottle
	Sensor error alarm
	Power failure alarm
	Door ajar alarm
	Low battery alarm
	Operation indicator
	Alarm validation (audible, visual, and remote alarm) Alarm mute with ring back
	Remote alarm contacts
	Four adjustable shelves (adjustable in ½" increments)
	Forced draft air circulation
	Keyed door lock
	Magnetic door gasket for positive seal
	Heavy duty door handle
	LED interior lights are safety shielded and switch controlled
	Swivel casters (3 ½")
	Probe access port (3/8")
	White powder-coated steel exterior and interior
	Three year parts and labor warranty, plus an additional two year compressor parts
	warranty
	Exterior dimensions: 28 3/8" W x 36 ½" D x 81 ¾" H with casters; 78 ¼" H without
	casters
	Cycle defrost
	115V, 60 Hz, 3 Amps, 1/5 HP
	ENERGY STAR CERTIFIED UL listed
	Shipping weight: 275 lbs.
	Onipping weight. 270 ibs.