Product Name:	ABT-HC-72S-TS
Description:	TempLog Premier Solid Door Laboratory Refrigerator
	72 Cu. Ft. capacity
	Three swing solid doors, self-closing, right hinged, not reversible
	Microprocessor temperature controller Adjustable operating temperature range: 1°C - 10°C
	TempData-2000 Touch Screen Digital Temperature Display and Alarm Module
	with battery back-up to support temperature monitoring and alarms
	Password-protected control panel
	Sample simulated digital temperature display; .1° Resolution, adjustable
	temperature display in Celsius or Fahrenheit
	Active temperature scale bar display
	Resettable 1-point min/max history
	Digital 24-hour temperature chart
	Temperature data logging
	USB port for data transfer and software updates
	Real-time operational display clock
	4-20mA output
	Electronic troubleshooting guide
	Electronic QR Code to access Owner's Manual & technical support
	contact information
	Adjustable alarm volume Audible and visual high/low temperature alarms
	2 probes; 1 in air and 1 in sample bottle
	Sensor error alarm
	Power failure alarm
	Door ajar alarm
	Low battery alarm
	Operation indicator
	Alarm validation (audible, visual, and remote alarm)
	Alarm mute with ring back
	Remote alarm contacts
	Twelve adjustable shelves (adjustable in ½" increments)
	Forced draft air circulation
	Keyed door locks
	Magnetic door gasket for positive seal Heavy duty door handles
	LED interior lights are safety shielded and switch controlled
	Swivel casters (3 ½")
	Probe access port (3/8")
	White powder-coated steel exterior and interior
	Hydrocarbon, natural refrigerant (R290)
	Three year parts and labor warranty, plus an additional two year compressor
	parts warranty
	Exterior dimensions: 81" W x 34" D x 81 3/4" H with casters; 78 1/4" H without
	casters
	Cycle defrost
	115V, 60 Hz, 7.5 Amps, 1/2 HP
	ETL listed
	Shipping weight: 600 lbs.