

Product Data Sheet

14 Cu. Ft. Premier Series Manual Defrost Freezer

ABT-HC-MFP-14P

ımage





Product Description

These manual defrost laboratory freezers offer professional cold storage for scientific and industrial settings. Units come equipped with temperature alarms, keyed door locks, probe access, and precise temperature monitoring through microprocessor control.

Our Premier Series freezers include premium additions such as extensive alarm systems and digital touch pad displays. All units run on natural, hydrocarbon refrigerant for environmental health and energy efficiency.

Certifications





General Description and Application Storage capacity (cu. ft) 14 Door Solid, right hinged, magnetic door gasket for positive seal 6 Fixed Shelves Freezer Compartments 7 Inner Doors Baskets Non-applicable Mounting Leveling Legs Interior lighting Non-applicable Airflow management Non-applicable External probe access Rear wall port (3/4") dia. Insulation Cabinet is foamed-in-place with EPA compliant high density urethane foam Exterior materials White powder coated steel Access Control Keyed door lock General warranty Two (2) year parts and labor warranty Compressor Parts Warranty Five (5) year compressor parts warranty Products require 4" of clearance around back and sides of the unit for proper ventilation - improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure Warranty Disclaimer Product Weight (lbs) 199

Shipping Weight (lbs) 219
Rated Amperage 3.5 Amps

Power Plug/Power Cord NEMA 5-15 Plug, conforms to UL471 requirements

Facility Electrical Requirement 115V, 60 Hz

Agency Listing and Certification ETL Listed, Energy Star Certifed

Refrigeration System

Refrigerant R290 Defrost Manual

Controller, Configuration, Alarms and Monitoring

Controller technology Microprocessor temperature controller with digital temperature display, °C/°F switchable

Display technology Bi-Colored LED display and alarm module with .1°C resolution, touchpad controls, resettable 1-point min/max history

Adjustable Temperature Range -15°C to -25°C

External alarm connection State switching remote alarm contacts

Alarms Audible and visual high and low temperature alarms, sensor error alarm, power failure alarm, alarm validation (audible, visual, and remote alarm), and alarm mute with ring back


Controller probe 2 Probes; 1 in air and 1 in sample bottle


Disclaimers Defrosting of the unit will need to be performed periodically to prevent severe ice build-up

Indoor use only, +18°C to +26°C (+65°F to +78°F), <70% RH

Specifications are subject to change

Dimensions						
	Width (in.)	Depth (in.)	Height (in.)	Door Swing (in.)	Total open Depth (in.)	
Exterior	23 7/8"	24 1/8"	74 3/4"	22 1/4"	46 1/4"	
Interior	16 1/2"	16 1/2"	Non-applicable			





Contact		
Customer Service	800-648-4041 Option 3	<u>customerservice@horizonscientific.com</u>
Technical Service	800-648-4041 Option 5, Parts Option 4	<u>technicalservice@horizonscientific.com</u>
Rev_05082023		