ABT-HC-MFP-17
Standard Natural Refrigerant Manual Defrost Laboratory
Freezer 17 Cu. Ft. Capacity
One solid swing door, right hinged
Microprocessor temperature controller
Digital temperature display
Audible and visual high and low temperature alarms
Remote alarm contacts
2 probes; 1 in air and 1 in glycol bottle
Adjustable operating temperature range: -15°C – -25°C
Six fixed shelves
Keyed door lock
Magnetic door gasket for positive seal
Leveling legs
Probe access port (3/4")
Powder coat exterior
High density urethane foam cabinet and door insulation
Natural, Hydrocarbon Refrigerant (R290)
Two year parts and labor warranty, plus additional three
year compressor parts warranty Exterior dimensions: 30 3/4" W x 28 5/8" D x 67 3/4"H
Note: This unit must have 4" clearance on sides and back for
adequate ventilation
Manual defrost
115V, 60 Hz, 5.5 Amps, ½ HP
ENERGY STAR CERTIFIED
ETL listed
Shipping weight: 220 lbs.