Product Name:	ABT-MFS-15-C
Description:	General Purpose Manual Defrost Laboratory Chest Freezer 15 Cu. Ft. capacity One solid door Mechanical thermostat Adjustable operating temperature range: -15°C – -25°C LED Interior light Two removable baskets Swivel casters (2 ¾") Keyed door lock Magnetic door gasket for positive seal White cabinet and lid High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging One-year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 56"W x 29 ½"D x 35 1/8"H with casters; 32 3/8" without casters (Height with hinged door open 90° is 62 ½") Manual defrost 115V, 60 Hz, 5 Amps, 1/4 HP Shipping weight: 197 lbs.