Product Name:	PH-ABT-HC-UCBI-0420-ADA-LH
Description:	PH-ABT-HC-UCBI-0420-ADA-LH Premier Pharmacy/Vaccine Undercounter Freezer Built-In-LEFT HINGED 4.2 Cu. Ft. Capacity One solid swing door, left hinged Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in sample bottle Adjustable operating temperature range: 5°F – -13°F (-15°C – -25°C) Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, "buffered" probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm Two fixed shelves
	Two fixed shelves
	Vaccine storage power cord warning label Keyed door lock
	Magnetic door gasket for positive seal
	Recessed door handle Leveling legs
	Probe access port (3/8")
	Powder coat exterior High density urethane foam cabinet and door insulation
	HFC-free refrigerant (R600a), foam insulation and packaging
	Pharmacy refrigerator/freezer toolkit and temperature logs
	Two year parts and labor warranty, plus an additional three year
	compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible
	Exterior dimensions: 23 3/4"W x 24 1/2"D x 31 15/16"H
	Manual defrost
	115V, 60 Hz, 1.5 Amps, 1/8 HP Energy Consumption 0.97 kWh/day
	ADA compliant
	UL/C-UL listed
	ENERGY STAR CERTIFIED
	Shipping weight: 133 lbs.