

Product Name:	PH-ABT-HC-UCBI-0420-ADA
Description:	<p>Premier Pharmacy/Vaccine Undercounter Freezer Built-In  4.2 Cu. Ft. Capacity  One solid swing door, right hinged*  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  Adjustable operating temperature range: 5°F – -13°F (-15°C – -25°C)  Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual &amp; audible temp alarm  Two fixed shelves  Vaccine storage power cord warning label  Keyed door lock  Magnetic door gasket for positive seal  Recessed door handle  Leveling legs  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 23 ¾ "W x 24 ½" D x 31 15/16 "H  Manual defrost  115V, 60 Hz, 1.5 Amps, 1/8 HP  Energy Consumption 0.97 kWh/day  ADA compliant  UL/C-UL listed  ENERGY STAR CERTIFIED  Shipping weight: 133 lbs.  *This product is available left hinged, product number PH-ABT-HC-UCBI-0420-ADA-LH</p>