

|               |   |
|---------------|---|
| Product Name: | PH-ABT-HC-UCFS-0430   |
| Description:  | <p>Premier Undercounter Pharmacy Freezer Freestanding (-30°C)<br/> 4 Cu. Ft. Capacity<br/> One solid swing door, right hinged, not reversible<br/> Microprocessor temperature controller<br/> Digital temperature display<br/> Audible and visual high and low temperature alarms<br/> 2 probes; 1 in air and 1 in sample bottle<br/> Remote alarm contacts<br/> Adjustable operating temperature range:-4°F to -22°F (-20°C to -30°C)<br/> Factory Preset at -30°C<br/> Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall<br/> Bottom footed shelf to allow air circulation<br/> Back wall evaporator cover<br/> Keyed door lock<br/> Magnetic door gasket for positive seal<br/> Recessed door handle<br/> Leveling legs<br/> LED interior light<br/> Probe access port (3/8")<br/> White cabinet and door<br/> High density urethane foam cabinet and door insulation<br/> HFC-free refrigerant (R290), foam insulation and packaging<br/> Two year parts and labor warranty, plus an additional three year compressor parts only warranty<br/> Pyxis®, compressor parts only Omnicell® and warranty AcuDose RX® compatible<br/> Exterior dimensions: 23 3/4" W x 24"D x 32 1/8"H<br/> <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i><br/> Manual defrost<br/> 115V, 60 Hz, 2.2 Amps, 1/4 HP<br/> ENERGY STAR CERTIFIED<br/> UL/C-UL listed<br/> Shipping weight: 128 lbs.</p> |