

Product Name:	PH-ABT-HC-UCFS-0440
Description:	<p>Premier Undercounter Pharmacy Freezer Freestanding (-40°C)  4 Cu. Ft. Capacity  One solid swing door, right hinged, not reversible  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature range: -35°F - -45°F (-37°C – -43°C)  <b>Please Note:</b> The unit temperature can be reset to -20°C  Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall  Bottom footed shelf to allow air circulation  Back wall evaporator cover  Keyed door lock  Magnetic door gasket for positive seal  Recessed door handle  Leveling legs  LED interior light  Probe access port (3/8")  White cabinet and door  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R290), foam insulation and packaging  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 23 ¾" W x 24"D x 32 1/8"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost  115V, 60 Hz, 2.2 Amps, 1/4 HP  ENERGY STAR CERTIFIED  UL/C-UL listed  Shipping weight: 128 lbs.</p>