

Product Name:	AJ-IV CS
Description:	<p>ABS Auto Jr. with Fiberboard Package System</p> <p>Total Vial Capacity (2 mL vials): 6000</p> <p>Total Box Capacity: 60</p> <p>Rack Capacity: 6</p> <p>Boxes per rack: 10</p> <p>Vial per box (2 mL): 81/100</p> <p>LN2 Capacity (Liters): 165</p> <p>Static Evaporation Rate (Liter/Day): 0.84</p> <p>Static Holding Time (days): 194</p> <p>Power requirement: 24 Vac</p> <p>Neck Diameter (mm): 215.9</p> <p>Exterior Height (mm): 990.6</p> <p>Exterior Diameter (mm): 683.3</p> <p>Sentry Controller</p> <p>Controller features:</p> <ul style="list-style-type: none"> ○ Temperature display (Celsius) ○ Level accuracy with 4 direct read sensors ○ Large format LED display ○ Alarm ring back ○ Redundant level sensor system ○ Remote monitoring system available ○ Liquid level indicator ○ Alarm package (audible, visual, and remote) including level (high and low), high temperature, and sensor failure <p>UL-Listed transformer</p> <p>Stainless Steel armor-plated transfer hose</p> <p>State of the art hanging rack design</p> <p>Cast aluminum roller base</p> <p>Lid with locking provision & protective boot</p> <p>Fiberglass restricted opening neck design</p> <p>One-year warranty on all parts and labor</p> <p>Three-year vacuum warranty</p> <p>Package System Includes:</p> <ul style="list-style-type: none"> - (60) 81 or 100 Cell Fiberboard Boxes & Dividers - (6) 10-2 Racks - (1) Pair CryoGloves - (1) LN₂ Measuring Ruler