Product Name:	AJ-IV
Description:	ABS Auto Jr.
	Total Vial Capacity (2 mL vials): 6000
	Total Box Capacity: 60
	Rack Capacity: 6
	Boxes per rack: 10
	Vial per box (2 mL): 81/100
	LN2 Capacity (Liters): 165
	Static Evaporation Rate (Liter/Day): 0.84
	Static Holding Time (days): 194
	Power requirement: 24 Vac
	Neck Diameter (mm): 215.9
	Exterior Height (mm): 990.6
	Exterior Diameter (mm): 683.3
	Sentry Controller
	Controller features:
	Temperature display (Celsius)
	Level accuracy with 4 direct read sensors
	Large format LED display
	Alarm ring back     Bedundent level concer system
	<ul> <li>Redundant level sensor system</li> <li>Remote monitoring system available</li> </ul>
	Remote monitoring system available     Liquid level indicator
	<ul> <li>Alarm package (audible, visual, and remote) including level</li> </ul>
	(high and low), high temperature, and sensor failure
	UL-Listed transformer
	Stainless Steel armor-plated transfer hose
	State of the art hanging rack design
	Cast aluminum roller base
	Lid with locking provision & protective boot
	Fiberglass restricted opening neck design
	One-year warranty on all parts and labor
	Three-year vacuum warranty