

Product Name:	AMIV-1 CS
Description:	<p>ABS Auto Max System with Fiberboard Package System  Auto Fill Tank  Total Vial Capacity (2 mL vials): 10,400  Total Bag Capacity: 235  LN2 Capacity (Liters): 165  Static Evaporation Rate (Liters/Day): 5  Normal working time (days): 33  Power requirement: 24 Vac  Exterior Width (mm): 586.7  Overall Height (mm): 1117.6  Exterior Depth (mm): 774.7  Usable Height (mm): 736.6  Internal Diameter (mm): 533.4  Sentry Controller</p> <p><b>Controller features:</b></p> <ul style="list-style-type: none"> <li>○ Temperature display (Celsius)</li> <li>○ Level accuracy with 4 direct read sensors</li> <li>○ Large format LED display</li> <li>○ Alarm ring back</li> <li>○ Redundant level sensor system</li> <li>○ Remote monitoring system available</li> <li>○ Liquid level indicator</li> <li>○ Alarm package (audible, visual, and remote) including level (high and low), high temperature, and sensor failure</li> </ul> <p>UL-Listed transformer  Stainless steel vacuum tank  Recessed casters  Fiberglass top  Full access opening  Leveling platform  Stainless steel armor-plated transfer hose  7-inch foam insulated hinged lid with locking provision, protective boot, and pneumatic lifters  One-year warranty on all parts and labor  Three-year vacuum warranty</p> <p><b>Package System Includes:</b></p> <ul style="list-style-type: none"> <li>- (7) 13-2 Tower Racks (81 or 100 cell boxes)</li> <li>- (4) 13-2 Mini Tower Racks (25 cell mini boxes)</li> <li>- (91) 81 or 100 Cell Fiberboard Boxes &amp; Dividers</li> <li>- (54) 25 Cell Mini Fiberboard Boxes &amp; Dividers</li> <li>- (1) LN<sub>2</sub> Measuring Ruler</li> <li>- (1) Pair CryoGloves</li> </ul>