

Product Name:	AMIV-2 CS
Description:	<p>ABS Auto Max System with Fiberboard Package System Auto Fill Tank Total Vial Capacity (2 mL vials): 24,050 Total Bag Capacity: 535 LN2 Capacity (Liters): 365 Static Evaporation Rate (Liters/Day): 7 Normal working time (days): 52 Power requirement: 24 Vac Exterior Width (mm): 863.6 Overall Height (mm): 1117.6 Exterior Depth (mm): 977.9 Usable Height (mm): 736.6 Internal Diameter (mm): 787.4 Sentry Controller</p> <p>Controller features:</p> <ul style="list-style-type: none"> ○ Temperature display (Celsius) ○ Level accuracy with 4 direct read sensors ○ Large format LED display ○ Alarm ring back ○ Redundant level sensor system ○ Remote monitoring system available ○ Liquid level indicator ○ Alarm package (audible, visual, and remote) including level (high and low), high temperature, and sensor failure <p>UL-Listed transformer Stainless steel vacuum tank Recessed casters Full access opening Fiberglass top Stainless steel armor-plated transfer hose Leveling platform 7-inch foam insulated hinged lid with locking provision, protective boot, and pneumatic lifters One-year warranty on all parts and labor Three-year vacuum warranty</p> <p>Package System Includes:</p> <ul style="list-style-type: none"> - (17) 13-2 Tower Racks (81 or 100 cell boxes) - (6) 13-2 Mini Tower Racks (25 cell mini boxes) - (221) 81 or 100 cell fiberboard boxes & dividers - (78) 25 cell mini fiberboard boxes & dividers - (1) LN₂ Measuring Ruler - (1) Pair CryoGloves