Product Name:	AMIV-2
Description:	ABS Auto Max System
	Auto Fill Tank
	Total Vial Capacity (2 mL vials): 24,050
	Total Bag Capacity: 535
	LN2 Capacity (Liters): 365
	Static Evaporation Rate (Liters/Day): 7
	Normal working time (days): 52
	Power requirement: 24 Vac
	Exterior Width (mm): 863.6
	Overall Height (mm): 1117.6
	Exterior Depth (mm): 977.9
	Usable Height (mm): 736.6
	Internal Diameter (mm): 787.4
	Sentry Controller
	Controller features:
	 Temperature display (Celsius)
	 Level accuracy with 4 direct read sensors
	 Large format LED display
	 Alarm ring back
	 Redundant level sensor system
	 Remote monitoring system available
	 Liquid level indicator
	 Alarm package (audible, visual, and remote) including level
	(high and low), high temperature, and sensor failure
	UL-Listed transformer
	Stainless steel vacuum tank
	Recessed casters
	Full access opening
	Fiberglass top
	Stainless steel armor-plated transfer hose
	Leveling platform
	7-inch foam insulated hinged lid with locking provision, protective
	boot, and pneumatic lifters
	One-year warranty on all parts and labor
	Three-year vacuum warranty