

|               |   |
|---------------|---|
| Product Name: | CRT-ABT-HC-S16S   |
| Description:  | <p>Upright Controlled Room Temperature Cabinet , Solid Door<br/> 16 Cu. Ft. capacity<br/> One swing solid door, self-closing, right hinged, not reversible<br/> Microprocessor temperature controller<br/> Digital temperature display<br/> Audible and visual high and low temperature alarms<br/> Remote alarm contacts<br/> Adjustable operating temperature range: 68°F–77°F (20°C – 25°C)<br/> Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual &amp; audible temp alarm<br/> 2 probes; 1 in air and 1 in sample bottle<br/> Three adjustable shelves<br/> Vaccine storage power cord warning label<br/> Forced draft air circulation<br/> Heavy duty door handle<br/> Keyed door lock<br/> Magnetic door gasket for positive seal<br/> LED interior lights are safety shielded and switch controlled<br/> Swivel Casters (3”)<br/> Probe access port (3/8”)<br/> White powder-coated steel exterior and interior<br/> Hydrocarbon, natural refrigerant (R290)<br/> One year parts and labor warranty, plus an additional four year compressor parts warranty<br/> Exterior dimensions: 25” W x 29 3/4” D x 79” H with casters; 75 1/2” H without casters<br/> Cycle defrost<br/> 115V, 60 Hz, 3 Amps, 1/5 HP<br/> ENERGY STAR CERTIFIED<br/> ETL listed<br/> Shipping weight: 277 lbs.</p> |