

Product Name:	ABT-FRP-04
Description:	<p>Standard Undercounter Flammable Storage Refrigerator  Freestanding 5 Cu. Ft. capacity  One solid swing door, right hinged, not reversible  Microprocessor temperature controller  Digital temperature display with battery back-up to support temperature monitoring and alarms  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature range: 1°-10°C  Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall  Keyed door lock  Magnetic door gasket for positive seal  Integrated door handle  One probe access port (3/8")  White cabinet and door  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Two year parts and labor warranty  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Exterior dimensions: 23 5/8" W x 24 ¼ " D x 32 3/16" H  Cycle defrost  115V, 60 Hz, 1.2 Amps, 1/9 HP  UL/C-UL listed  Units meet NFPA guidelines 45 and 70  Shipping weight: 118 lbs.</p>