

Product Name:	PH-ABT-HC-SSP-49FA
Description:	<p>Premier Series Pharmacy Stainless Steel Freezer 49 Cu. Ft. capacity  Two swing solid doors, self-closing, not reversible  6 ea. Adjustable Shelves  Digital display microprocessor temperature controller  Adjustable operating temperature: 5°F – -13°F (-15°C - -25°C)  Temperature Monitor Device - complies with the current CDC guidelines, 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory, °F/°C switchable, field installable, and visual &amp; audible temperature alarms  2 probes; 1 in air and 1 in sample bottle  Probe attached to evaporator for defrost termination sensing  Sensor error alarm  Continuous product (simulated) temperature display  Remote alarm contacts  Forced draft air circulation  Keyed door lock  Magnetic door gaskets for positive seal  Integrated door handle  LED interior lights are safety shielded and switch controlled  Probe access port (3/4")  Stainless steel exterior  White-powder coated aluminum interior with 304 stainless steel floor  Easily cleaned smooth interior  Hydrocarbon natural refrigerant (R290)  2 year parts and labor warranty, plus an additional 3 year compressor parts warranty  Auto defrost  4 swivel casters (5 3/4")  115V/60 Hz  ETL listed  Exterior dimensions: 54” W x 35 3/4” D x 83 1/4” H with casters; 77 1/2” H without casters  115V, 60 Hz, 9 Amps, 3/4 HP  ETL listed  Shipping weight: 495 lbs.</p>