

Product Name:	PH-ABT-HC-SSP-49
Description:	<p>Premier Pharmacy/Vaccine Stainless Steel Refrigerator  49 Cu. Ft. Capacity  Two stainless steel swing doors  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  2 probes; 1 in air and 1 in sample bottle  Adjustable operating temperature range: 36°F – 46°F (2°C- 8°C)  Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Six adjustable shelves  Vaccine storage power cord warning label  Forced draft air circulation  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Swivel casters (adds 6 ½" to height)  Interior light  Probe access port (3/8")  Stainless steel interior floor  Stainless steel exterior  High density urethane foam cabinet insulation  Hydrocarbon, natural refrigerant (R290)  Pharmacy refrigerator/freezer toolkit and temperature logs  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 54"W x 35"D x 84"H with casters; 77 ½"H without casters  Cycle defrost  115V, 60 Hz, 4.5 Amps, 1/4 HP  ETL listed  Shipping weight: 450 lbs.</p>