

Product Name:	ABT-HC-RFC9-LH
Description:	<p>Premier Refrigerator & Freezer Combination-LEFT HINGED 9 Cu. Ft. capacity (4.6 Ref./4.2 Freezer) Two exterior solid swing doors, left hinged Microprocessor temperature controller, refrigerator and freezer Digital temperature display, refrigerator and freezer Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C Audible and visual high and low temperature alarms, refrigerator and freezer 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer Remote alarm contacts, refrigerator and freezer 2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall Back wall evaporator cover, Refrigerator Forced draft circulation Refrigerator Keyed door locks Magnetic door gasket for positive seal Recessed door handle Leveling legs Stacking kit Two probe access ports (3/8"); one, refrigerator/one, freezer White powder coated exterior HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 23 ¾"W x 24 ½"D x 66 ¾"H Manual defrost freezer / cycle defrost refrigerator 115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer UL/C-UL listed Shipping weight: 269 lbs. *Please note this unit requires 2 outlets.</p>