

Product Name:	ABT-HC-UCBI-0404-ADA-LH
Description:	<p>Premier Undercounter Refrigerator Built-In 4.6 Cu. Ft. Capacity One solid swing door, left hinged Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in glycol bottle Adjustable operating temperature range: 1°C – 10°C Three shelves (two adjustable/one fixed) with guard rail on back to prevent product from pushing against back wall Back wall evaporator cover Forced draft circulation Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis[®], Omnicell[®] and AcuDose RX[®] compatible</p> <p>Exterior dimensions: 23 ¾ "W x 24 ½" D x 31 15/16 "H Cycle defrost 115V, 60 Hz, 1.74 Amps, 1/8 HP ENERGY STAR CERTIFIED ADA Compliant and UL/C-UL Listed Shipping weight: 150 lbs.</p>