

Product Name:	PH-ABT-HC-SSP-49
Description:	Premier Pharmacy/Vaccine Stainless Steel Refrigerator 49 Cu. Ft. Capacity Two stainless steel swing doors Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in sample bottle Adjustable operating temperature range: 36°F – 46°F (2°C- 8°C) Calibrated using a NIST traceable device, certificate included Six adjustable shelves Vaccine storage power cord warning label Forced draft air circulation Keyed door locks Magnetic door gasket for positive seal Integrated door handles Swivel casters (adds 6 ½" to height) Interior light Probe access port (3/8") Stainless steel interior floor Stainless steel exterior High density urethane foam cabinet insulation Hydrocarbon, natural refrigerant (R290) Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 54"W x 35"D x 84"H with casters; 77 ½"H without casters Cycle defrost 115V, 60 Hz, 4.5 Amps, 1/4 HP ETL listed Shipping weight: 450 lbs.