

Product Name:	PH-ABT-HC-UCBI-0420
Description:	Premier Pharmacy/Vaccine Undercounter Freezer Built-In 4.2 Cu. Ft. Capacity One solid swing door, right hinged* Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in glycol bottle Adjustable operating temperature range: 5°F – -13°F (-15°C – -25°C) Calibrated using a NIST traceable device, certificate included Two fixed shelves Vaccine storage power cord warning label Keyed door lock Recessed door handle Magnetic door gasket for positive seal Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 23 ¾ "W x 24 ½" D x 33 3/8 "H Manual defrost 115V, 60 Hz, 1.5 Amps 1/8 HP Energy Consumption 0.97 kWh/day UL/C-UL listed ENERGY STAR CERTIFIED Shipping weight: 132 lbs. *This product is available left hinged, product number PH-ABT-HC-UCBI-0420-LH