

Product Name:	PH-ABT-HC-UCFS-0440
Description:	Premier Undercounter Pharmacy Freezer Freestanding (-40°C) 4 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: -35°F - -45°F (-37°C – -43°C) Please Note: The unit temperature can be reset to -20°C Calibrated using a NIST traceable device, certificate included Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall Bottom footed shelf to allow air circulation Back wall evaporator cover Keyed door lock Magnetic door gasket for positive seal Recessed door handle Leveling legs LED interior light Probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R290), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 23 ¾" W x 24"D x 32 1/8"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost 115V, 60 Hz, 2.2 Amps, 1/4 HP ENERGY STAR CERTIFIED UL/C-UL listed Shipping weight: 128 lbs.