

<p>Product Name:</p>	<p>ABT-FFP-30-TS</p>
<p>Description:</p>	<p>TempLog Premier Flammable Storage Freezer  30 Cu. Ft. capacity  One solid swing door, right hinged, not reversible  Microprocessor temperature controller  Adjustable operating temperature range: -15°- -25°C  Touch Screen Digital Temperature Display and Alarm Module  Password-protected control panel  Sample simulated digital temperature display; .1° Resolution, temperature display in Fahrenheit  Active temperature scale bar display  Resettable 1-point min/max history  Digital 24-hour temperature chart  Temperature data logging  USB port for data transfer and software updates  Real-time operational display clock  4-20mA output  Electronic troubleshooting guide  Electronic QR Code to access Owner's Manual &amp; technical support contact Information  Adjustable alarm volume  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Door ajar alarm  Low battery alarm  Operation indicator  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  Remote alarm contacts  Four fixed shelves  Keyed door lock  Magnetic door gasket for positive seal  Stainless steel door handle  Leveling legs  One probe access port (3/8")  White textured cabinet and door  High density urethane foam cabinet and door insulation  CFC-free refrigerant (R404a), foam insulation and packaging  Three year parts and labor warranty  Exterior dimensions: 35 3/4" W x 35" D x 76 3/4"H  Manual defrost  115V, 60 Hz, 6 Amps, 1/3 HP  UL/C-UL listed  Units meet NFPA guidelines 45 and 70  Shipping weight: 418 lbs.</p>