

Product Name:	ABS 4
Description:	<p> ABS CryoMax System Total Vial Capacity: 6000 (2ml) Box Capacity: 60 Rack Capacity: 6 Vials per box: 81/100 Liquid Nitrogen Capacity: 175 liters Static Evaporation Rate: 0.99 liters per day Static Holding Time: 104 days Neck Diameter (mm): 215.9 Exterior Height (mm): 958.8 Exterior Diameter (mm): 660.4 Durable aluminum tank with powder coated exterior Stainless steel storage racks, configured for 2" high boxes Lid lock provision Insulated neck tube core Compatible design for both liquid and vapor phase storage One-year warranty on all parts and labor Three-year vacuum warranty </p>