

Product Name:	ABS LD 10
Description:	ABS LD Class 10 Liter Liquid Nitrogen Capacity 0.18 Liter/Day Static Evaporation Rate Neck Diameter (mm): 56 Exterior Height (mm): 546.1 Exterior Diameter: 261.6 Designed for efficient storage of LN2 Low LN2 consumption Durable composite lid High strength neck tube reduces LN2 loss Strong, lightweight aluminum construction Safer operation & maintenance One-year parts & labor warranty Three-year vacuum warranty