


Mid Arm Length CryoGloves

ABS CG MA S WP, ABS CG MA M WP, ABS CG MA L WP, ABS CG MA XL WP
ABS CG MA S, ABS CG MA M, ABS CG MA L, ABS CG MA XL

CUTSHEET



Size	A	B	C	D
XL/11	10"	5.5"	6"	7"
L/10	9.5"	5.5"	5.5"	6.5"
M/9	9"	5"	5"	6"
S/8	8.5"	5"	4.5"	5.5"

Sizing: CryoGuard gloves are available in two styles (waterproof and standard water repellent) and four sizes: Small, Medium, Large and Extra Large. These four sizes correspond to European and American sizes of 8, 9, 10 and 11. Sizing tolerance of gloves are $\pm 1/4"$.

Disclaimer: The gloves must be loose so they can be thrown off quickly if liquid spills into them. Under no circumstances should any of these gloves be submerged directly in liquid nitrogen or used with liquid oxygen.